

## Post Processing

Equipment	Chemicals
slide rack and dishes (Shandon Lipshaw #121)	1XSSC
500 ml beaker, 2 stir bars and shaker	1M Sodium Borate, filtered (make from Boric acid, adjust pH with NaOH)
hot plate and 4L beaker	1-Methyl-2-pyrrolidinone (anhydrous, 99.5%, FW=99.13)
humid chamber (Sigma H6644) and heating block	Succinic Anhydride (99+%, FW=100.07, a moisture sensitive irritant)
centrifuge with adapter for slide rack.	95% Ethanol

### Before you start

- a) heat up heat block to at least 90°C and invert the insert
  - b) place slide-staining chambers in front of the heating block.
  - c) fill slide staining chamber with SSC or water.
  - d) get a dust-free board needed for the cross-linking step.
  - e) place in the chem. hood:
    - large glass beaker with about 1 L boiling water, cover with foil as it is heated
    - stirring plate for preparation of blocking solution
    - shaker
  - f) for the preparation of the blocking solution, have a clean 500 ml beaker ready.
  - g) check if the following chemicals are available:
    - succinic anhydride (6.0 g)
    - 1-Methyl-2-pyrrolidinone (350 ml)
    - 15 ml 1 M Na-Borate, pH 8.0 (filtered)
1. Rehydrate the slides
    - a) Rehydrate slides by inverting (array side down toward steamy water) them over warm water in a slide-staining chamber. Let the spots become glistening and juicy.  
*Be careful not to allow the water to touch the array.*
    - b) Immediately flip them (array side up!) onto a heating block (inverted about 90°C). Watch the steam evaporate. When the array spots dry in a rapid wave-like pattern, remove them from the heating block. This takes about 5 seconds. Do 1 slide at a time.
  2. UV cross link
    - a) Place the slides, array side up, on a flat, dust-free board that fits into the UV cross-linker.  
*Do not put them on a saran wrap surface since the slides stick to it.*
    - b) Irradiate with 600 uJoules UV light.

3. At this time, we place the slides in the slide rack and secure it with a strip of metal wire on top. Otherwise, there is a good chance that the slides will pop out of the rack during the vigorous washing in step 7.

*Skip prewash steps #4-#6 if the spots on the array are small. This step is to reduce pluming for arrays with an excess amount of dna.*

- 4. Gently wash slides with 1x SSC / 0.01%SDS for 30 sec**
- 5. Gently wash slides with 0.06x SSC for 30 sec.**
- 6. Dry slides in centrifuge at 600rpm 1-2 min.**

7. Block free lysine (wear a lab coat when working with methyl pyrrolidinone)
  - a) Add 6.0 g of succinic anhydride into 350 ml 1-Methyl-2-pyrrolidinone while stirring.
  - b) As soon as the solids dissolve, quickly add 15 ml of the 1M Na-Borate pH 8, and pour the mixed solution into a slide washing tray (make sure the tray is small enough that 350ml of solution will completely cover the slides when they are placed inside).
  - c) Quickly place the slides into the succinic anhydride solution and vigorously plunge for 60 seconds. Rotate at 60 rpm for 15 min.

NOTE: The slides should be rapidly plunged into the organic acetic anhydride solution to quickly wash the unattached DNA away from the slides.

**Do not pour the solution over the slides.**

- d) Remove the slide rack from the organic reaction mixture and place it immediately into the boiling water bath for 90 seconds.
- e) Transfer the slide rack to the 95% ethanol wash tray and carry it directly to the tabletop centrifuge.
- f) Spin dry the slides by centrifugation at 600 rpm for 2 min. Use a counter balance with the same number of slides in a rack.
- g) Carefully transfer the slides to a dry slide box for storage in a dessicator.